




Helping Families, One Bed At A Time



FEATURES

WHAT MAKES THIS BED DIFFERENT...

- PATENTED SECURE CLAMPING DESIGN FOR ENCLOSURE PANELS (#7,434,280)
- PATENTED SECURE ZIPPER WITH ENCLOSURE (#10,959,495)
- INDEPENDENT ENCLOSURE PANELS
- TENSION CONTROL TABS
- LIFT ACCESSIBILITY
- EASY TRANSFER HEIGHT FROM MATTRESS
- 8-INCH THICK MATTRESS
- 3 YEAR WARRANTY



SPECIFICATIONS

STANDARD WITH EVERY BED...

BEDFRAME:

- SOLID HARDWOOD – NO METAL
- SOLID MATTRESS PLATFORM
- NO ARTICULATION
- SLEEPING AREA DIMENSIONS:
80”L x 36”W x 51”H (Standard)
80”L x 60”W x 51”H (Queen)
- 3 STAIN COLORS – NO UPCHARGE

ENCLOSURE PANELS:

- HEAVY-WEIGHT CANVAS
- WATER RESISTANT/REPELLANT
- WINDOWS – HEAVY-WEIGHT MESH
- INDEPENDENT PANELS
- PADDED VERTICAL CORNERS
- OPEN TOP – PANEL AVAILABLE
- 12 COLORS – NO UPCHARGE

MATTRESS:

- PRESSURE REDISTRIBUTION FOAM
- REMOVABLE, FLUID RESISTANT & FIRE RETARDANT COVER
- TOP OF MATTRESS HEIGHT FROM FLOOR IS 24”

ENCLOSED SAFETY BED

A SOLUTION FOR PROTECTION

FDA REGISTERED & LISTED

AUTHORIZED COVERAGE UNDER;

MEDICAID, PRIVATE INSURANCES,
WAIVER PROGRAMS

ADDITIONAL OPTIONS

CAN BE ADDED TO ANY BED

TOP MESH PANEL

- THIS OPTION WILL FULLY ENCLOSE BED
- CAMERA ACCESS CUT-OUTS CAN BE ADDED

TUBING PORTS

- PROVIDE ACCESS INTO SLEEPING AREA WITHOUT COMPROMISING THE SAFETY AND SECURITY OF ENCLOSURE
- INSTALLED IN HEAD AND FOOT PANELS

EXTERIOR PADDING

- PROVIDES PADDING TO OUTSIDE BED FRAME
- 12 COLORS – NO UPCHARGE

ZIPPERED BACK PANEL

- PROVIDES A SECOND ACCESS INTO SLEEPING AREA

LOW BED HEIGHT

- PROVIDES OPTION TO CRAWL INTO BED
- LOWERS MATTRESS HEIGHT TO 13"
- SLEEPING AREA REMAIN SAME SIZE

RETRACTABLE CASTERS

- PROVIDES ABILITY TO EASILY MOVE BED

QUEEN SIZE

- CONVERSION KIT AVAILABLE

COLOR OPTIONS

ENCLOSURE PANELS & EXTEIOR PADDING



Aruba



Black



Burgundy



Capri



Charcoal




Erin Green



Logo Red



Mediterranean
Blue




Pink



Sea



Toast



Natural

WE UNDERSTAND THE IMPORTANCE OF CHOOSING THE RIGHT BED.

PLEASE FEEL FREE TO REACH OUT WITH ANY QUESTIONS.

COLOR OPTIONS

BED FRAME STAIN

NATURAL ● HONEY ● WALNUT



MADE IN MAINE WITH 100% USA SOURCED MATERIALS



WOMEN OWNED SMALL BUSINESS CERTIFIED