

# DEBUG ALUMINUM ALL-TERRAIN WALKER



# DEBUG ALUMINUM ALL-TERRAIN WALKER



The Debug Aluminum All-Terrain Walker provides stability to individuals walking on sandy beaches, rough terrain and even snow that normal walkers cannot accommodate. Light weight, portable and ready to go when you are.

LIGHT WEIGHT ALUMINUM FRAME FOLDS FOR TRANSPORTATION AND STORAGE / HAND BRAKE AND LOCKING MECHANISM / FOLDING SEAT AND BACK-PAD TO SIT / INCLUDED STORAGE BASKETS / STAINLESS STEEL CASTERS AND AXLES RESIST SALT AND SAND / PNEUMATIC WHEELS ROLL ACROSS SAND, SNOW AND GRAVEL / FITS THROUGH A STANDARD ADA DOORWAY / LOW MAINTENANCE DELRIN BUSHINGS / ADJUSTABLE AND REMOVABLE PUSH HANDLES

## SPECIFICATIONS

SPEC	ENGLISH (INCH)	METRIC (CM)
Overall height (adjustable)	34-39	86-99
Floor to seat height	26	66
Width	23	58
Chairweight	18 lbs	8 kg
Capacity	200 – 300 lbs	91 – 136 kg

## ADDITIONAL NOTES

---

Frame material: anodized aluminum

---

Walker frames are available in pediatric, adult and large adult sizes

---

Castors and axles are made from #316L stainless steel with 10 year warranty

---

Rollator folding style walker is able to fit through a standard ADA doorway and be used indoors and outside in areas where using a standard walker is difficult.

Manufactured in the USA,