



Home Bed Side Helper

Item #	Product Description	UOM
15065R-1	Retail	1/cs

- Attaches to Hollywood-style home bed frame with four legs on the floor. (Figure A)
- Designed to provide assistance getting in and out of home bed.
- Easy to install, no tools required.
- Rail features four easily accessible grasping points.
- Allows for normal bedding drape. (Figure B)
- Adjustable leg height to fit most beds. 11" wide frame for one-or two-handed use.
- Durable, solid steel tubing.
- Comes with 8" mounting bracket for a secure fit.
- Retail Packaged.

specifications

Rail Dimensions	11" (W) x 37" (H) x 11.5" (D)
Weight Limit	250 lbs.
Warranty	Limited Lifetime