



Full-Electric Bariatric Bed, 48"

Item #	Product Description	UOM
15302	Bariatric Bed	1/cs
15302BV-1HR	Bariatric Bed with 1 pair "T" Rails (item #15021BV)	1/cs
15302BV-PKG	Bariatric Bed Package (contains Bariatric Bed, 48" Bariatric Foam Mattress, item #15310, and 1 pair "T" Rails, item #15021BV)	1/cs

- Quiet, smooth operation.
- Heavy-duty frame ensures strength and patient safety.
- All steel construction.
- Larger sleep surface than conventional bed.
- Split-pan design with removable bed ends is easy to set up.
- No tools necessary for set up.
- Hand control (included) provides for multiple bed positioning for patients.
- Emergency manual crank included.
- Low air loss/alternating pressure mattress available.
- "T" rails available. (Sold separately) (Figure A)
- Weight Capacity 750 lbs.

specifications

Bed Height	18.25" - 26.25"
Unit Dimensions	48" (W) x 88" (L)
Weight	312 lbs.
Carton Shipping Weight	343 lbs.
Weight Limit	750 lbs.