



FIGURE A

- Delivers a uniform oxygen pulse at increments: 1, 2, 3, 4, 5 and 6
- Two continuous flow (CF) settings preset @ 2 LPM and 4 LPM
- Delivers up to a 6:1 savings ratio from 14-40 breaths per minute
- Senses a breath/delivers oxygen in the first half of the inspiratory cycle
- Operating pressure range: 200-3000 PSI
- Light weight and quiet operation
- Easy-to-operate single selector knob
- Compact design fits easily in virtually any oxygen carry bag

PRODUCTS**Description**

OM-812	Bonsai Velocity Conservator, 1/ea
--------	-----------------------------------

SPECIFICATIONS**Description**

Warranty	Limited 2 Year
----------	----------------

ACCESSORIES**Description**

SEAL100	Seal Washer,Conservators/Regulators,10/pkg, PKG
---------	---