



COMPARE TO		
drive™	BREWER®	CLINTON®
13045	11100	M22



Mayo-Instrument Stand

Double Post

Item #	UOM
13045	1/cs

- Double post provides added support for use in surgical procedures.
- Removable stainless steel tray measures 19" x 12.5".
- Adjustable height from 32.5"– 52".
- Two hooded casters.

specifications	
Tray Dimensions	19" (W) x 12.5" (D)
Height	32.5" - 52"
Base Dimensions	21" (W) x 19" (D)
Weight	11 lbs
Carton Shipping Weight	14 lbs.
Warranty	Limited Lifetime

accessories		
13035T	Replacement Tray	1/cs

3-Panel Privacy Screen

Item #	UOM
13508	1/cs

- Sturdy, 1" diameter anodized aluminum tubing.
- Screens are 6 mil. flame retardant bacteriostatic white vinyl.
- Four, 3" hooded casters.
- Easy to maintain and clean.
- Each hinged panel measures: 56" (H), 27 ½" (W)

specifications	
Overall Height	70 ½"
Weight	15 lbs.
Carton Shipping Weight	20 lbs.
Warranty	Limited Lifetime