



### 3-Position Recliner, Bariatric Extra Wide

Item #	Color		UOM
D574EW-BR	Blue Ridge		1/cs
D574EW-J	Jade		1/cs
D574EW-R	Rosewood		1/cs

- 3 positions: upright (Figure A), deep recline and elevated leg rest/footrest. (Figure B)
- Retractable lock bar secures the chair in the desired position
- Self-adjusting suspension in the seat prevents bottoming out.
- Side panels “pop off” for easy cleaning and maintenance.
- Large, blow-molded tray locks in 5 positions and can be stored on the side of the chair.
- 5" casters (2 with locks) make it easy to move and secure the recliner.
- Wide, comfortable armrests are angled to provide support in any reclined position.
- Neutral-colored, sturdy steel frame is easy to clean and maintain.
- Moisture barrier on seat prevents seepage.
- Comfortable built-in headrest.
- Foam and vinyl meet California Flammability Standard and Boston Fire Department specifications.

specifications	
Overall Dimensions	47.5" (H) x 32" (W) x 41.5" (D)
Width Between Arms	24.5"
Arm Height From Floor	26.5"
Arm Height From Seat	9.5"
Seat Dimensions	21" (H) x 25" (W) x 20" (D)
Carton Shipping Weight	127 lbs.
Weight Limit	500 lbs.
Warranty Limited Lifetime on Frame, 1 Year on Upholstery and Moving Parts	