



Free-Standing Trapeze

with Base (Comparable to Invacare®'s Bariatric Trapeze, model #BARTRAP)

Item #	Product Description	UOM
13049SV	For use by patients up to 650 lbs	1/cs

- Designed to assist larger individuals in changing positions while in bed and aid in transferring in and out of bed.
- Easy-to-assemble tool-free design.
- Lightweight design for easy delivery.
- Chain is adjustable from edge of headboard to the end of boom by sliding to desired position.

specifications

Total Dimensions	72" (H) x 31" (W) x 51" (D)
Inside Width Between Legs	26"
Outside Width Between Legs	31"
Boom Length	37"
Weight	74 lbs.
Carton Shipping Weight	83 lbs.
User Weight Limit	650 lbs.

Free-Standing Silver Vein Trapeze

with Base and Wheels

Item #	Product Description	UOM
13039SV	For use by patients up to 1000 lbs	1/cs

- Designed to assist larger individuals to change positions while in bed and aid in the transfer from bed to chairs.
- Easy to assemble. (Requires only a drill, rubber mallet and socket with wrench)
- The overhead boom swivels from 0° to 180°.
- Chain is adjustable from edge of headboard to the end of boom by sliding to desired position.

specifications

Total Dimensions	73" (H) x 43" (W) x 48.5" (D)
Inside Width Between Legs	37.5"
Outside Width Between Legs	48"
Boom Length	36"
Weight	100 lbs.
Carton Shipping Weight	110 lbs.
User Weight Limit	1000 lbs.