



Overdoor Exercise Pulley

Cervical Traction Set

Item #	UOM
13005	1/bx

- Designed to safely and easily increase range of motion.
- Easily fits over standard door.
- Comes complete with door bracket, pulleys, cord and handles.
- For use when sitting or standing.
- Retail packaged.

specifications	
Warranty	Limited Lifetime

accessories		
13005P	Pulley Only	1/ea
13005R	Rope Only	1/ea
13005H	Handle Only	1/ea

Item #	UOM
13004	1/bx

- Heavy-duty head halter comes complete with metal support and self-attaching closures.
- Complete with 12" adjustable spreader bar, 8' traction rope, double-sealed rings, water bag and "S" hook.
- Machine washable.
- One size fits all.
- Retail packaged.

specifications	
Warranty	Limited Lifetime

accessories		
13004H	Head Halter Only	1/ea
13004W	Water Bag Only	1/ea
13005P	Pulley Only	1/ea
13004R	Rope Only	1/ea