



(Shown with Fusion Maple Top Finish and U Base)



Non-Tilt Premium Multi-Feature Overbed Table

- Pneumatic cylinder lift mechanism provides infinite adjustment from 27" - 42".
- High-pressure laminate top with raised urethane lip for spill containment.
- Dual removable cup holders with flexible grip tabs.
- Storage between cup holders – perfect for pens, cell phones, etc.
- Dual hand grips and 2" casters for easy maneuvering.
- Column can be mounted on left or right side.
- Available with "H" (Figure A) or "U" base in chrome or opal powder coat.
- Optional dual-mirror vanity drawer.
- Base and columns available in additional powder-coat finishes.
- Available on drop ship basis only.



16005 H

TOP FINISH	H Base Chrome Frame		H Base Opal Frame		UOM
	W/O VANITY	W/VANITY	W/O VANITY	W/VANITY	
Gunstock Walnut	16005-HCGW	16005-HCGWV	16005-HOGW	16005-HOGWV	1/cs
Fusion Maple	16005-HCFM	16005-HCFMV	16005-HOFM	16005-HOFMV	1/cs
Banister Oak	16005-HCBO	16005-HCBOV	16005-HOBO	16005-HOBOV	1/cs
Wild Cherry	16005-HCWC	16005-HCWCV	16005-HOWC	16005-HOWCV	1/cs
English Oak	16005-HCEO	16005-HCEOV	16005-HOEO	16005-HOEOV	1/cs
Figured Mahogany	16005-HCDM	16005-HCDMV	16005-HODM	16005-HODMV	1/cs
Nepal Teak	16005-HCNT	16005-HCNTV	16005-HONT	16005-HONTV	1/cs
Quartered Oak	16005-HCQO	16005-HCQOV	16005-HOQO	16005-HOQOV	1/cs
Solar Oak	16005-HCSO	16005-HCSOV	16005-HOSO	16005-HOSOV	1/cs

16005U

TOP FINISH	U Base Chrome Frame		U Base Opal Frame		UOM
	W/O VANITY	W/VANITY	W/O VANITY	W/VANITY	
Gunstock Walnut	16005-UCGW	16005-UCGWV	16005-UOGW	16005-UOGWV	1/cs
Fusion Maple	16005-UCFM	16005-UCFMV	16005-UOFM	16005-UOFMV	1/cs
Banister Oak	16005-UCBO	16005-UCBOV	16005-UOBO	16005-UOBOV	1/cs
Wild Cherry	16005-UCWC	16005-UCWCV	16005-UOWC	16005-UOWCV	1/cs
English Oak	16005-UCCEO	16005-UCCEOV	16005-UOEO	16005-UOEOV	1/cs
Figured Mahogany	16005-UCDM	16005-UCDMV	16005-UODM	16005-UODMV	1/cs
Nepal Teak	16005-UCNT	16005-UCNTV	16005-UONT	16005-UONTV	1/cs
Quartered Oak	16005-UCQO	16005-UCQOV	16005-UOQO	16005-UOQOV	1/cs
Solar Oak	16005-UCSO	16005-UCSOV	16005-UOSO	16005-UOSOV	1/cs

specifications

Height (adjustable)	27" - 42"
Weight	51 lbs.
Carton Shipping Weight	55 lbs.
Weight Limit	100 lbs.
Base Dimensions	18" x 32"
Tabletop Dimensions	18" x 32"
Warranty	Limited Lifetime