



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Tool-free assembly of seat and legs (Figure A)
- Drainage holes in bench reduce slipping
- Aluminum frame is lightweight, durable and corrosion proof
- Legs are height adjustable in 1" increments, and angled with suction-style tips to provide additional stability
- Made using recycled materials

PRODUCTS	Description
RTL12203KDR	Bath Seat w/o Back, KD, Retail, 1 EA
SPECIFICATIONS	Description
Carton Shipping Weight	7 lbs
Outside Legs	19.5" (W) x 19.5" (D)
Seat Dimensions	19.25" (W) x 11.5" (D)
Seat Height	14 21
Warranty	Limited Lifetime
Weight	5 lbs 5 lbs 4.5 lbs
Weight Capacity	300 lbs