



Portable Shower Bench

Bath Mat

Item #	Color	UOM
12023	White	1/cs

Item #	Product Size	Color	UOM
12950	Large	White	1/bx

- One-piece shower bench adds convenience, comfort and safety to the bathroom.
- It is portable for easy storage and travel.
- Handle provides additional safety and stability for the user.
- Rubber stops help protect tub from scratches.
- Drain holes prevent pooling.

- Mat is held securely by multiple suction cups.
- Extra-long bath mat adds safety and security by providing a large, slip-proof surface in the bath.
- Soft rubber mat is easy to clean and roll up for storage.
- Mat has a contour-cut design to fit around drain.
- Mold resistant latex-free rubber.

specifications

Width	14.5"
Length	30.5"
Weight	5.5 lbs.
Carton Shipping Weight	7 lbs.
Weight Limit	250 lbs.

specifications

Width	15.75"
Length	35.5"
Weight	2.4 lbs.
Carton Shipping Weight	3.6 lbs.
Weight Limit	Unlimited