



A



B

Portable Bath Step

Item #	Product Description	UOM
RTL12068	Contains 1 Bath Step per box	1/cs

- Makes getting in and out of the bath tub easier and safer.
- Unique modular design allows the multiple steps to be stacked. (Figure A)
- Step height is 4". When stacking, each additional step adds 2".
- Bath step comes with four removable covers on the surface. When stacking, remove covers to allow steps to lock together. (Figure B)
- Attractive, textured surface is slip resistant and comfortable.
- Retail Packaged.

specifications

Dimensions	19.5" (W) x 4" (H) x 16" (D)
Weight Limit	450 lbs.