




12487 Shown



Folding Shower Chair

Item #	Product Description	UOM
12486	Without Back	1/cs
12487	With Back	1/cs

- Product available with back (main image) and without. (Figure A)
- Conveniently folds flat. (Figure B)
- Blow-molded bench provides comfort and strength.
- Drainage holes in seat reduce slipping.
- Aluminum frame is lightweight, durable and corrosion proof.
- Angled legs provide additional stability.
- Retail packaged.

specifications

Outside Legs (12486)	20.5" (W) x 19.25" (D)
Outside Legs (12487)	17.75" (W) x 14" (D)
Seat Dimensions (12486)	17.5" - 18.5"
Seat Dimensions (12487)	8" (W) x 11" (D)
Seat Height (12486 12487)	17.5" - 18.5" 16.25" - 18.25"
Weight (12486 12487)	4.5 lbs. 5 lbs.
Carton Shipping Weight (12486 12487)	6.6 lbs. 7 lbs.
Weight Limit	300 lbs.

Premium Series Shower Chair with Back and Arms

Item #	Product Description	UOM
RTL12505	With Back and Arms	4/cs

- Easy, tool-free height adjustment with height indicator. Adjusts in ½" increments. (Figure A)
- Enhanced comfort seat fits all shapes and sizes.
- Easy, tool-free removable back, legs and arms.
- Lightweight design has large stable legs with soft non-slip feet.
- Retail packaged.

specifications

Outside Legs	17" (W) x 15.5" (D)
Seat Dimensions	21.75" (W) x 16.5" (D)
Seat Height	16" (L) x 20.5" (H)
Weight	10 lbs.
Carton Shipping Weight	48 lbs.
Weight Limit	350 lbs.