



RTL10370RC
Red Crackle



Aluminum Folding Canes with Gel Grip, Height Adjustable

"T" Handle Cane, Designer Colors

Item #	Color	UOM
RTL10370RC	Red Crackle	1/cs
RTL10370BC	Blue Crackle	1/cs
RTL10370BK	Black	1/cs

Item #	Color	UOM
10335BKW-1	Black Wave	1/ea

- New gel grip reduces hand stress and fatigue. (Figure A)
- Aluminum folding cane, height adjustable.
- Cane folds into four convenient parts for easy storage.
- Manufactured with sturdy, extruded aluminum tubing.
- Comes with plastic clip to hang cane when folded.
- Available in Red Crackle (main image), Blue Crackle (Figure B) or Black. (Figure C)
- Handle height adjusts from 32" to 36".
- 300 lb. Weight Capacity.

- Attractive wood handle with brass collar.
- Comes standard with wrist strap.
- Manufactured with sturdy, extruded aluminum tubing.
- Easy-to-use, one-button height adjustment with locking ring. Locking ring prevents rattling. (Figure A)
- Handle height adjusts from 30" to 39".
- 300 lb. Weight Capacity.

specifications

Handle Height	32" - 36"
Weight Limit	300 lbs.
Warranty	Limited Lifetime

specifications

Handle Height	30" - 39"
Weight Limit	300 lbs.
Warranty	Limited Lifetime