




RTL10310
Black



Quad Canes, Small Base

Item #	Product Description	Color
10301-4	Vinyl Contoured Grip	Chrome
10301F-4	Foam Rubber Grip	Chrome
RTL10310	Foam Rubber Grip	Black
10311-4	Vinyl Contoured Grip	Bronze
10312FP-1	Foam Rubber Grip	Floral

- Four point base, combined with offset handle provides additional stability and support.
- Easy-to-use, one-button height adjustment with locking bolt, allows quick, safe fitting.
- Available in Black (main image), Chrome (Figure A), Bronze (Figure B) or Floral. (Figure C)
- No assembly required.
- Handle height adjusts from 30" to 39".
- Handle available in foam rubber (main image) or vinyl. (Figure D)
- 300 lb. Weight Capacity.

specifications

Handle Height	30" - 39"
Weight Limit	300 lbs.