



## Deluxe Steel Drop-Arm Commode Knocked-Down Frame

Item #	UOM
11125KD-1	1/cs
11125KD-2	2/cs

- Easy-to-assemble frame.
- Easy-to-release drop-arm mechanism allows for safe lateral patient transfers to and from commode.
- Durable plastic snap-on seat and lid.
- Easy-to-clean gray powder-coated steel finish.
- Plastic armrests and back provide extra comfort and support.
- Comes with 12 qt. commode bucket with carry handle and splash shield.

### specifications

Width	19.75" between arms   23.5" Outside
Seat Dimensions	13.75" (W) x 16.5" (D)
Height	20" - 24"
Weight	18 lbs.
Carton Shipping Weight (1   2)	20   39 lbs.
Weight Limit	300 lbs.

### accessories

RTL12085	Biodegradable Sanitary Commode Liner	12 liners/box, 6bx/cs
----------	--------------------------------------	-----------------------