



Never Clean a Pail Again!



A



Pre-Moistened Wash Glove

Biodegradable Sanitary Commode Liner

Item #

RTL12086

- New Wash Glove is the ideal solution for waterless bathing.
- Made of highly resistant, soft, non-woven embossed material.
- The wash glove is pre-moistened with a non-rinse, protective washing lotion without soap or alcohol, leaving the skin clean and dry.
- The wash glove is hypoallergenic and dermatologically tested.
- Safe to use on all types of skin from newborns to seniors.
- Comes in a resealable pack and can be microwaved for more comfort.
- Single use.
- Ready to use.

Item #

RTL12085

- The new Drive commode liner is a single-use medical device to be positioned over a commode.
- The liner is equipped with a super-absorbent pad which turns bodily waste (up to 17 oz.) into gel within seconds.
- The liner is easy to use, convenient and safe.
- Prevents the soiling/contamination of commode pails and reduces unpleasant odors.
- Safe and secure collection, transport and disposal of contaminated bodily fluids.
- Eliminates odors.
- Avoid spills and splashes.
- Time and money saving. No hand washing.
- Provides comfort and dignity to the patient and patient's family.