



A



Dressing Stick

Item #	Product Description	Product Size	UOM
RTL2030	Wood	18"	
RTL2032	Padded, with Shoe Horn	24"	

- Both 18" and 24" models have a plastic coated "S" hook at one end. RTL2030 has a small "C" hook on the other.
- Can be used for both pulling on and pushing off clothes. (Figure A)
- Retail packaged.

specifications

Warranty Limited 6 Months

Molded Stocking Aid

Item #	UOM
RTL2010	

- Makes it possible to put on stockings without bending.
- Curved cutouts on the sides hold the sock while it is being pulled up the leg.
- Stocking aid holds the sock open and allows the user to easily pull the sock or stocking up with minimal effort.
- Retail packaged.

specifications

Warranty Limited 6 Months