



EZ Assist Pole with Rotating Handle

Padded Swivel Seat Cushion

Item

15500TH

- Assists an individual to and from a sitting position in a chair, on a bed or in the bathroom.
- Rotating handle pivots and locks in multiple positions and is height adjustable.
- Rotating handle is 18" in width.
- Multi-level rails allow for gripping in multiple positions.
- Sturdy, steel construction with brown-vein finish is rust resistant.
- 300 lb. Weight Capacity.

specifications

Fits Ceiling Heights	100.25" - 126.75"
Weight Capacity	300 lbs.
Floor to Bottom of Handle	Lowest: 25" Highest: 36.5"
Warranty	Limited Lifetime

Item

RTLAGE-300

- Swivel seat allows individual to turn up to 360 degrees.
- Makes getting on or off a seat or chair easy and comfortable. (Figure A)
- Durable and stable, with a 300 lb. Weight Capacity.
- Portable: weighs only 3.25 lbs.
- Dimensions: 15 5/8" diameter.
- Retail Packaged.

specifications

Weight Capacity	300 lbs.
Warranty	Limited 1 Year