

Bariatric Sentra EC Heavy-Duty, Extra-Extra-Wide Wheelchair
Dual Cross Brace



HCPCS
K0007

Shown with
optional PH-SF
swing-away footrests

Item #	Seat Width	Armrest	Product Front Riggings	UOM
STD26ECDDA	26"	Detachable Desk Arm	Sold Separately	1/cs
STD26ECDF A	26"	Detachable Full Arm	Sold Separately	1/cs
STD28ECDDA	28"	Detachable Desk Arm	Sold Separately	1/cs
STD28ECDF A	28"	Detachable Full Arm	Sold Separately	1/cs
STD30ECDDA	30"	Detachable Desk Arm	Sold Separately	1/cs
STD30ECDF A	30"	Detachable Full Arm	Sold Separately	1/cs

- 700 lb. Weight Capacity.
- Carbon steel frame with triple-coated chrome for an attractive, chip-proof, maintainable finish.
- Reinforced steel gussets at all weight bearing points provide additional strength.
- Durable reinforced nylon upholstery.
- Steel, spoke wheels are strong and maintenance free.
- Heavy-duty, flat-free 8" front casters are adjustable in two positions.
- Extra-heavy-duty front forks.
- Reinforced side frame, caster journals and dual cross bars.
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability.
- Extra-large, padded armrests provide comfort and support.
- Riggings have extra-large, heavy-duty, aluminum foot plates that provide maximum support, are attractive and crack proof. (Must be ordered separately)
- Comes with push-to-lock wheel locks.

specifications

Width Open Closed with 26" Seat	34.5" 14"	Product Weight with 20" Seat	62 lbs.
Width Open Closed with 28" Seat	37" 14.5"	Product Weight with 22" Seat	64 lbs.
Width Open Closed with 30" Seat	39" 15"	Product Weight with 24" Seat	66 lbs.
Height	35"	Carton Shipping Weight with 20" Seat	68 lbs.
Seat-To-Floor Height	19.5"	Carton Shipping Weight with 22" Seat	70 lbs.
Back Height	18"	Carton Shipping Weight with 24" Seat	72 lbs.
Seat Depth	20"	Weight Limit	700 lbs.
Length with Riggings	43"		

accessories

Front Riggings

PH-SF	Aluminum Foot plate	1/pr
PHELR	Aluminum Foot plate	1/pr

Anti-Tippers with Wheels

STDS829	Anti-Tippers With Wheels	1/pr
---------	--------------------------	------

Wheel Lock Extension

STDS801	Wheel Lock Extension	1/ea
---------	----------------------	------

Bariatric Seat Belts

STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric	1/ea
STDS856	Seat Belts, Velcro® Type Closure, Bariatric	1/ea

Telescoping IV Pole Attachment

STDS820	Telescoping IV Pole Attachment	1/ea
---------	--------------------------------	------

Universal Oxygen "E" Cylinder Holder

STDS803	Universal Oxygen "E" Cylinder Holder	1/ea
STDS804	Universal Oxygen "E" Cylinder Holder with IV Pole Attachment	1/ea

Wheelchair Cane/Crutch Holder

STDS1034	Wheelchair Cane/Crutch Holder	1/ea
----------	-------------------------------	------

Heel Loops

STDS831	Heel Loops	1/pr
---------	------------	------