



Weights only 19 lbs.



Deluxe Fly-Weight Aluminum Transport Chair with Removable Casters

Item #	Seat Width	Color		UOM
RTLFW19RW-BL	19"	Blue Frame, Black Upholstery		1/cs
RTLFW19RW-RD	19"	Red Frame, Black Upholstery		1/cs

- Weighs only 19 lbs.
- Aluminum frame for lightweight strength and durability.
- Easy folding and storage in closet or trunk. (Figure A)
- Composite wheels are lightweight and maintenance free.
- 8" wheels in front and rear easily remove without tools. (Figure B - C)
- Deluxe back release. (Figure D)
- Deluxe, all-aluminum rear wheel locks.
- Comes with carry pocket on backrest. (Figure E)
- Both 17" and 19" come with tool-free, swing-away footrests.
- Comes with seat belt for added safety.
- Flared arms allow the back to be folded virtually flat. (Figure F)
- Soft, plush nylon upholstery.
- Back folds down for easy storage and transport.
- Convenient cup holder can easily attach to chair. Sold separately. (Model# STDS1040S)
- Weight Capacity: 300 lbs.

specifications

Overall Width	22" (open) 9" (closed)
Overall Height	37.25"
Seat To-Floor Height	19"
Back Height	18"
Seat Dimensions	19" (W) x 16" (D)
Overall Length w/Riggings	33"
Weight	19 lbs.
Carton Shipping Weight	24 lbs.
Weight Limit	300 lbs.

accessories

STDS1040S	Universal Cup Holder	1/ea
-----------	----------------------	------