

## Lynx Ultra Lightweight Wheelchair



HCPCS: K516FBADDA-SF, K518FBADDA-SF, K520FBADDA-SF;  
K0005; E0951; E0973 | K516FBADDA-ELR,



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F




FIGURE G



FIGURE H

- Ultra Lightweight Aluminum Frame
- Flip-back, Padded, Height Adjustable Desk Arms Provides Added Patient Comfort
- Comes with push to lock wheel locks
- Treaded tires
- Quick release axles
- Dual axle positions
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability
- Adjustable seat depth
- Rear anti-tippers standard
- Swing away footrests (K516FBADDA-SF, K518FBADDA-SF, K520FBADDA-SF)
- Elevating Leg Rests (K516FBADDA-ELR, K518FBADDA-ELR, K520FBADDA-ELR)
- Mobility has never been easier with the new and improved lightweight Lynx K5 Wheelchair. Gone are the days of struggling with heavy and bulky wheelchairs. Instead, choose a wheelchair that keeps you moving and adjusts to your individual needs.

### PRODUCTS

	Description
K516FBADDA-ELR	K5 WC 16" F/B,ADJ,DeskArm, ELR
K516FBADDA-SF	K5 WC 16" F/B,ADJ,DeskArm, SF
K518FBADDA-ELR	K5 WC 18" F/B,ADJ,DeskArm, ELR
K518FBADDA-SF	K5 WC 18" F/B,ADJ,DeskArm, SF
K520FBADDA-ELR	K5 WC 20" F/B,ADJ,DeskArm, ELR
K520FBADDA-SF	K5 WC 20" F/B,ADJ,DeskArm, SF

### SPECIFICATIONS

	Description
Anti-Tip Wheels	Yes
Back of Chair Height	16.5"
Brakes	Push to lock
Casters	7"
Color	Black
Folded Dimensions	12.2" x 35.5" x 43.9"
Handle Height Range	35.5" - 37.5"
Height Setting (armrest to floor)	28.5"-31.1"
Number of Wheels	4
Rear Wheels	24"
Seat Depth	16"-18"
Seat To Floor Height	19.5"-17.5"

**SPECIFICATIONS****Description**

Seat Width	K516FBADDA-SF and K516FBADDA-ELR: 16"   K518FBADDA-SF and K518FBADDA-ELR: 18"   K520FBADDA-SF and K520FBADDA-ELR: 20"
Warranty	Limited Lifetime
Weight Capacity	275 lbs
Width Closed	12.2"
Material Composition	Aluminum; Nylon