


**Rebel Wheelchair**  
Single Axle



Item #	Seat Width	Armrest	Product Front Riggings	UOM
RTLREB18DDA-SF	18"	Flip-Back Desk Arm	Swing-away Footrest	1/cs

- Quick-release rear wheels allow for easy storage or transport. (Figure A)
- Fold-down back and folding frame allows transport in almost any trunk. (Figure B)
- Easily removable swing-away footrests are also convenient for storage.
- Tool-free removable front casters.
- Desk-style armrests make transfers easy and accessibility better.
- Lightweight and extremely maneuverable, easy to navigate tight areas.
- Large accommodating seating surface.
- Carton size makes the Rebel perfect for retail.



specifications			
Seat Width	18"	Overall Length (with riggings)	43"
Overall Width	25.5" (open)   10.5" (closed)	Weight	37 lbs.
Overall Height	35"	Carton Shipping Weight	42 lbs.
Seat-To-Floor Height	19.5"	Weight Limit	300 lbs.
Back Height	17"		
Seat Dimensions	18" (W) x 16" (D)		

accessories			
<b>Seat Belts</b>			
STDS850	Seat Belts, Auto-Clasp Type Closure	1/ea	
STDS851	Seat Belts, Velcro® Type Closure	1/ea	
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric	1/ea	
STDS856	Seat Belts, Velcro® Type Closure, Bariatric	1/ea	
<b>Universal Oxygen "E" Cylinder Holder</b>			
STDS803	Universal Oxygen "E" Cylinder Holder	1/ea	
STDS804	Universal Oxygen "E" Cylinder Holder with IV Pole Attachment	1/ea	
<b>Wheel Lock Extension</b>			
STDS801GT	Wheel Lock Extension	1/ea	
<b>Anti-Fold Bar</b>			
STDS806	Anti-Fold Bar	1/ea	
<b>Wheelchair Cane/Crutch Holder</b>			
STDS1034	Wheelchair Cane/Crutch Holder	1/ea	
<b>Heel Loops</b>			
STDS831	Heel Loops	1/pr	
<b>Cup Holder</b>			
STDS1040S	Universal Cup Holder	1/ea	