

## Silver Sport 2 Wheelchair


Dual axle



HCPCS  
K0001 / K0002

Item #	Color	Seat Width	Armrest	Product Front Riggings	UOM
SSP216FA-SF	Silver Vein	16"	Fixed Arm	Swing-away Footrests	1/cs
SSP216FA-ELR	Silver Vein	16"	Fixed Arm	Swing-away Elevating Leg rests	1/cs
SSP216DDA-SF	Silver Vein	16"	Detachable Desk Arm	Swing-away Footrests	1/cs
SSP216DDA-ELR	Silver Vein	16"	Detachable Desk Arm	Swing-away Elevating Leg rests	1/cs
SSP216DFA-SF	Silver Vein	16"	Detachable Full Arm	Swing-away Footrests	1/cs
SSP216DFA-ELR	Silver Vein	16"	Detachable Full Arm	Swing-away Elevating Leg rests	1/cs
SSP218FA-SF	Silver Vein	18"	Fixed Arm	Swing-away Footrests	1/cs
SSP218FA-ELR	Silver Vein	18"	Fixed Arm	Swing-away Elevating Leg rests	1/cs
SSP218DDA-SF	Silver Vein	18"	Detachable Desk Arm	Swing-away Footrests	1/cs
SSP218DDA-ELR	Silver Vein	18"	Detachable Desk Arm	Swing-away Elevating Leg rests	1/cs
SSP218DFA-SF	Silver Vein	18"	Detachable Full Arm	Swing-away Footrests	1/cs
SSP218DFA-ELR	Silver Vein	18"	Detachable Full Arm	Swing-away Elevating Leg rests	1/cs
SSP220DDA-SF	Silver Vein	20"	Detachable Desk Arm	Swing-away Footrests	1/cs
SSP220DDA-ELR	Silver Vein	20"	Detachable Desk Arm	Swing-away Elevating Leg rests	1/cs
SSP220DFA-SF	Silver Vein	20"	Detachable Full Arm	Swing-away Footrests	1/cs
SSP220DFA-ELR	Silver Vein	20"	Detachable Full Arm	Swing-away Elevating Leg rests	1/cs

- Powder-coated silver vein steel frame. (Figure A)
- Embossed vinyl upholstery is durable, lightweight, attractive and easy to clean.
- Urethane tires mounted on composite wheels provide durability and low maintenance.
- 8" front casters, adjustable in three positions.
- Padded armrests provide additional comfort.
- Carry pocket on backrest provides additional convenience.
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings.
- Plastic foot plates with heel loops on footrests.
- Comes with push-to-lock wheel locks.



### specification

Width Open with 16" Seat	24"	Product Weight with 18" Seat	41 lbs.
Width Open with 18" Seat	26"	Product Weight with 20" Seat	42 lbs.
Width Open with 20" Seat	28"	Carton Shipping Weight with 16" Seat	48 lbs.
Width Closed	12.5"	Carton Shipping Weight with 18" Seat	49 lbs.
Height	36"	Carton Shipping Weight with 20" Seat	50 lbs.
Seat-To-Floor Height	17.5" - 19.5"	Weight Limit with 16" Seat	250 lbs.
Back Height	16"	Weight Limit with 18" Seat	300 lbs.
Seat Depth	16"	Weight Limit with 20" Seat	350 lbs.
Length with Riggings	42"		
Product Weight with 16" Seat	40 lbs.		

### accessories

#### Anti-Tippers

STDS833N	Anti-Tippers With Wheels	1/pr
STDS818	Anti-Tippers Without Wheels	1/pr

#### Seat Belts

STDS850	Seat Belts, Auto-Clasp Type Closure	1/ea
STDS851	Seat Belts, Velcro® Type Closure	1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric	1/ea
STDS856	Seat Belts, Velcro® Type Closure, Bariatric	1/ea

#### Telescoping IV Pole Attachment

STDS820	Telescoping IV Pole Attachment	1/ea
---------	--------------------------------	------

#### Universal Oxygen "E" Cylinder Holder

STDS803	Universal Oxygen "E" Cylinder Holder	1/ea
STDS804	Universal Oxygen "E" Cylinder Holder with IV Pole Attachment	1/ea

#### Wheel Lock Extension

STDS801GT	Wheel Lock Extension	1/ea
-----------	----------------------	------

#### Overhead Anti-Theft Device

STDS823	Overhead Anti-Theft Device	1/ea
STDS821	Overhead Anti-Theft Device With IV Hooks	1/ea

#### Anti-Fold Bar

STDS806	Anti-Fold Bar	1/ea
---------	---------------	------

#### Wheelchair Cane/Crutch Holder

STDS1034	Wheelchair Cane/Crutch Holder	1/ea
----------	-------------------------------	------

#### Heel Loops

STDS831	Heel Loops	1/pr
---------	------------	------

#### Limb Support

WASR	Limb Support, Right	1/ea
WASL	Limb Support, Left	1/ea

#### Chart Carry Pocket

STDS836	Chart Carry Pocket	1/ea
---------	--------------------	------

#### Universal Forward Stabilizer

STDS862	Universal Forward Stabilizer	1/pr
---------	------------------------------	------