



## Steel Transport Chair

Item #	Seat Width	Color	UOM
TR37E-SV	17"	Silver Vein Frame, Black Upholstery	1/cs
TR39E-SV	19"	Silver Vein Frame, Black Upholstery	1/cs

- Weighs 26 lbs.
- Attractive silver vein finish. (Figure A)
- Available in 17" or 19" seat widths.
- Back folds down for storage and transport. (Figure B)
- Swing-away footrests come standard.

- Padded armrests provide additional comfort.
- 8" casters in front and rear.
- Comes with attractive nylon upholstery that is durable and easy to clean.
- Standard rear wheel locks.
- Comes with upholstery handles to fold chair.
- Convenient cup holder can easily attach to chair. Sold separately. (Model# STDS1040S)

### specifications

Overall Width (TR37)	21.5" (open)   8" (closed)
Overall Width (TR39)	23.5" (open)   10.75" (closed)
Overall Height	39"
Seat-To-Floor Height	19"
Back Height	18"
Seat Dimensions (TR37)	17" (W) x 16" (D)
Seat Dimensions (TR39)	19" (W) x 16" (D)
Overall Length w/Riggings	33"
Weight	26 lbs.
Carton Shipping Weight	31 lbs.
Weight Limit	250 lbs.

### accessories

STDS1040S	Universal Cup Holder	1/ea
-----------	----------------------	------