



Weighs only 14.5 lbs.



## TranSport Aluminum Transport Chair

Item #	Seat Width	Color	UOM
TS19	19"	Blue Upholstery, Blue Frame <span style="color: blue;">●</span>	1/cs

- Weighs 14.5 lbs.
- Aluminum frame is lightweight and strong.
- Frame magnets keep footrests secure against frame to prevent swinging open when the chair is carried.
- Padded flip-back arms are removable.
- Upholstery is comfortable, lightweight and easy to clean.
- Two-way back release allows to push or pull lever.
- The transport chair conveniently folds for travel and storage. (Figure A)
- Back folds forward.
- Comes with tool-free height adjustable swing-away footrests.
- Aluminum rear wheel locks.
- Weight Capacity: 250 lbs.

specifications	
Overall Width	24.5" (W) (open)   10" (W) (closed)
Overall Height	36"
Seat-To-Floor Height	18"
Back Height	16"
Seat Dimensions	19" (W) x 16" (D)
Overall Length w/ Riggings	38"
Weight	14.5 lbs.
Carton Shipping Weight	19.5 lbs.
Weight Limit	250 lbs.