

## Silver Sport Recliner with Vinyl Upholstery

HCPCS:



- Vinyl upholstery is comfortable, lightweight and easy-to-clean
- Hydraulic reclining mechanism allows for infinite adjustments up to 180°
- Arms have support extension for reclining position
- Carbon Steel frame with silver vein finish is attractive, and easy-to-maintain
- Composite Mag-style wheels are lightweight and maintenance free
- 8" front casters
- Padded armrests provide added patient comfort
- Wheels set back on frame prevent tipping
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Rear anti-tippers standard
- Rear anti-tippers standard headrest extension with cushioned head immobilizer standard
- Comes standard with swing-away elevating leg rests
- Metal foot plates
- Comes with push to lock wheel locks

**PRODUCTS****Description**

SSP16RBDDAV	SilverSport16"ReclBack,DetDesk,Vnl,1/ea
SSP16RBDFAV	SilverSport16"ReclBack,FullArm,Vnl,1/ea
SSP18RBDDAV	SilverSport18"ReclBack,DetDesk,Vnl,1/ea
SSP18RBDFAV	SilverSport18"ReclBack,FullArm,Vnl,1/ea
SSP20RBDDAV	SilverSport20"ReclBack,DetDesk,Vnl,1/ea
SSP20RBDFAV	SilverSport20"ReclBack,FullArm,Vnl,1/ea

**SPECIFICATIONS****Description**

Back Height	21"   31"
Carton Shipping Weight	66 lbs   67 lbs   67.5 lbs
Height	50"
Length with Riggings	50"
Material	Steel; Vinyl; Plastic
Product Weight with 16 Seat	59 lbs
Product Weight with 20 Seat	61 lbs
Product Weight with 18 Seat	60 lbs
Seat Depth	16"
Seat To Floor Height	20.5"
Warranty	Limited Lifetime
Weight Limit with 16 Seat	300 LBS

**SPECIFICATIONS****Description**

Weight Limit with 18 Seat	300 LBS
Weight Limit with 20 Seat	350 LBS
Width Closed	12"
Width Open	23"   25"   28"
Arm Rest Type (s)	Detachable Desk Arm   Detachable Full Arm

**ACCESSORIES****Description**

STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDS1040S	Cup Holder, Swivel,1BX,6/CS, BX
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
STDS818	Anti-Tippers,rear,SentraRecliner,1pr/bx, BX
STDELR-AL	ArticulatingElev.Legrest, Std,1pr/cs, CS
LK3JELR	ELR Set,Black,1Pair eq EA, PR
WASR	Limb Support, Right,1/cs, CS