

NEW

\*\*\*\*\*

1P  
ACK



## Flex N Go Cane

Item #	Product Description	UOM
RTL10305	Flex N Go Cane	1/ea

- Tip flexes and swivels to maintain contact with the ground when the cane is at an angle, providing superior balance.
- Newly designed tip is removable without tools, for use with most other canes.
- Height adjusts for use by most individuals.
- 3 point base allows the cane to stand on its own.
- Conveniently folds for easy storage.

### specifications

Handle Height	32.5" – 39.5"
Number of Height Adjustments	8
Weight Limit	300 lbs.