

bath bench comes
completely assembled



WEIGHT CAPACITY: 300 LBS.

instructions

Periodically check assembly for continued tightness. **NEVER USE BATH SAFETY SEAT IF ANY LOOSENESS IN THE ASSEMBLY IS NOTICED, AS PERSONAL INJURY MAY RESULT. SIT. DO NOT STAND** on this seat to avoid a possible dangerous fall.

CAUTION: Maximum weight limitation for this product not to exceed 300 lbs.

Meets ANSI Standard A117.1 Specifications.

NOTE: To maintain the beauty of the original finish, clean with mild soap and rinse with warm water. Wipe bath seat dry after each use.