

GRAVIS Floor Lift, Manual Base 600 lbs

HCPCS: E0635



- Weight capacity of 600 lbs
- Base option: Manual
- Elegant, contemporary design with soft lines but still heavy duty
- Oversized sealed rear casters roll easily and lock securely
- Easily lifts resident from floor
- Redundant controls allow for use even if hand control is inoperable
- Smart & Safety Technology LCD display Visual battery fuel indicator on Hand Pendant and Lift Control Module Auditable low battery alert fuel indicator Lift safety stop if resident or caregiver accidentally gets in the way of the boom arm Audible Alert if the lift detects weight overload Maintenance monitoring: by months or lift cycles Audible Alert that maintenance is due
- Visual battery fuel indicator on Hand Pendant and Lift Control Module
- Auditable low battery alert fuel indicator
- Lift safety stop if resident or caregiver accidentally gets in the way of the boom arm
- Audible Alert if the lift detects weight overload
- Maintenance monitoring: by months or lift cycles
- Audible Alert that maintenance is due
- NEW Six-point cradle with spring clips offer a more secure and wider foot print reducing sling squeeze, shallow height
- offers greater height range
- Can be charged on, or off the lift
- Includes 1 battery and external charger stand

PRODUCTS	Description
FLNP600	Gravis Floor lift,manual base,600 lbs

SPECIFICATIONS	Description
Battery Pack Weight	7.9 lbs (5.0A), 4.9 lbs (2.9A)
Battery Warranty	1 yr
Front Caster Diameter	3"
Gross Weight	156.5 lbs
Ground Clearance	2.9"
Net Weight	132.2 lbs
Warranty on Battery	1 yr
Warranty on Electronics	5 yr

ACCESSORIES	Description
FLSCALE	Patient Lift Scale 660lbs capacity, EA
13012	Patient Sling,Canvas,1/bx, BX

ACCESSORIES

Description

13019	Chains for Patient Lft,1pr/bx, BX
13023	Patient Lift,New Style,1/cs, CS
13025	Sling,Polystr.for Lift 1/ea, BX
13026	Sling,Polystr.w/Cmde Open1/ea, BX
13057	Replacement Cradle for 13240, 13244 & FL, EA
13060	BariatricSling,Canvas,1/bx, BX
13061	BariatricSlingw/ComdeOpen,Canvas,1/bx, BX
13232	Straps for Patient Slings,1set, BX
13220L	U-Sling 4pt padded Lg.1/bx, BX
13220M	U Sling 4pt padded Med,1/bx, BX
13220S	U-Sling 4pt padded Sm,1/bx, BX
13220XL	U-Sling 4pt padded X-lg,1/bx, BX

GRAVIS Floor Lift, Power Base 600 lbs

HCPCS: E0635



- Weight capacity of 600 lbs
- Base option: Power
- Elegant, contemporary design with soft lines but still heavy duty
- Oversized sealed rear casters roll easily and lock securely
- Easily lifts resident from floor
- Redundant controls allow for use even if hand control is inoperable
- Smart & Safety Technology LCD display Visual battery fuel indicator on Hand Pendant and Lift Control Module Auditable low battery alert fuel indicator Lift safety stop if resident or caregiver accidentally gets in the way of the boom arm Audible Alert if the lift detects weight overload Maintenance monitoring: by months or lift cycles Audible Alert that maintenance is due
- Visual battery fuel indicator on Hand Pendant and Lift Control Module
- Auditable low battery alert fuel indicator
- Lift safety stop if resident or caregiver accidentally gets in the way of the boom arm
- Audible Alert if the lift detects weight overload
- Maintenance monitoring: by months or lift cycles
- Audible Alert that maintenance is due
- NEW Six-point cradle with spring clips offer a more secure and wider foot print reducing sling squeeze, shallow height offers greater height range
- Can be charged on, or off the lift
- Includes 1 battery and external charger stand

PRODUCTS	Description
FLP600	Gravis Floor lift, power base, 600 lbs

SPECIFICATIONS	Description
Battery Pack Weight	7.9 lbs (5.0A), 4.9 lbs (2.9A)
Front Caster Diameter	3"
Gross Weight	156.5 lbs
Ground Clearance	2.9"
Net Weight	132.2 lbs
Warranty on Battery	1 yr
Warranty on Electronics	5 yr

ACCESSORIES	Description
FLSCALE	Patient Lift Scale 660lbs capacity, EA
13012	Patient Sling,Canvas,1/bx, BX
13019	Chains for Patient Lft,1pr/bx, BX

ACCESSORIES

Description

13023	Patient Lift,New Style,1/cs, CS
13025	Sling,Polystr.for Lift 1/ea, BX
13026	Sling,Polystr.w/Cmde Open1/ea, BX
13057	Replacement Cradle for 13240, 13244 & FL, EA
13060	BariatricSling,Canvas,1/bx, BX
13061	BariatricSlingw/ComdeOpen,Canvas,1/bx, BX
13232	Straps for Patient Slings,1set, BX
13220L	U-Sling 4pt padded Lg.1/bx, BX
13220M	U Sling 4pt padded Med,1/bx, BX
13220S	U-Sling 4pt padded Sm,1/bx, BX
13220XL	U-Sling 4pt padded X-lg,1/bx, BX