

LEVANTAR Floor Lift, Manual Base 500 lbs

HCPCS: E0635



- Weight capacity of 500 lbs
- Base option: Manual
- Elegant, contemporary design with soft lines but still heavy duty
- Oversized sealed rear casters roll easily and lock securely
- Easily lifts resident from floor
- Redundant controls allow for use even if hand control is inoperable
- Smart & Safety Technology LCD display Visual battery fuel indicator on Hand Pendant and Lift Control Module Auditable low battery alert fuel indicator Lift safety stop if resident or caregiver accidentally gets in the way of the boom arm Audible Alert if the lift detects weight overload Maintenance monitoring: by months or lift cycles Audible Alert that maintenance is due
- Visual battery fuel indicator on Hand Pendant and Lift Control Module
- Auditable low battery alert fuel indicator
- Lift safety stop if resident or caregiver accidentally gets in the way of the boom arm
- Audible Alert if the lift detects weight overload
- Maintenance monitoring: by months or lift cycles
- Audible Alert that maintenance is due
- NEW Six-point cradle with spring clips offer a more secure and wider foot print reducing sling squeeze, shallow height offers greater height range
- Can be charged on, or off the lift
- Includes 1 battery and external charger stand

PRODUCTS

Description

FLNP500	Levantar Floor lift,manual base,500 lbs
---------	---

SPECIFICATIONS

Description

Battery Pack Weight	7.9 lbs (5.0A), 4.9 lbs (2.9A)
Front Caster Diameter	3"
Gross Weight	156.5 lbs
Ground Clearance	2.9"
Net Weight	132.2 lbs
Warranty on Battery	1 yr
Warranty on Electronics	5 yr

ACCESSORIES

Description

FLSCALE	Patient Lift Scale 660lbs capacity, EA
13012	Patient Sling,Canvas,1/bx, BX
13019	Chains for Patient Lft,1pr/bx, BX

ACCESSORIES

Description

13025	Sling,Polystr.for Lift 1/ea, BX
13026	Sling,Polystr.w/Cmde Open1/ea, BX
13060	BariatricSling,Canvas,1/bx, BX
13061	BariatricSlingw/ComdeOpen,Canvas,1/bx, BX
13232	Straps for Patient Slings,1set, BX
13220L	U-Sling 4pt padded Lg.1/bx, BX
13220M	U Sling 4pt padded Med,1/bx, BX
13220S	U-Sling 4pt padded Sm,1/bx, BX
13220XL	U-Sling 4pt padded X-lg,1/bx, BX
13221L	Sling w/Cmde.Mesh,Lg,1/bx, BX
13221M	Sling,w/Comde,MeshMed,1/bx, BX
13221XL	Sling w/Cmde,Mesh,XL,1/bx, BX
13221XXL	Sling Full Body XXL Mesh 1/bx, BX
13222L	Sling Full Body Lg Solid 1/bx, BX
13222M	Sling Full Body Med Solid 1/bx, BX
13223L	Sling Full Body Lg Mesh 1/bx, BX
13223M	Sling Full Body Med Mesh 1/bx, BX

LEVANTAR Floor Lift, Power Base 500 lbs

HCPCS: E0635



FIGURE A




FIGURE B



FIGURE C




FIGURE D



FIGURE E




FIGURE F



FIGURE G



FIGURE H

- Weight capacity of 500 lbs
- Base option: Power
- Elegant, contemporary design with soft lines but still heavy duty
- Features you expect in a Premium lift
- Oversized sealed rear casters roll easily and lock securely
- Easily lifts resident from floor
- Redundant controls allow for use even if hand control is inoperable
- Smart & Safety Technology LCD display Visual battery fuel indicator on Hand Pendant and Lift Control Module Auditable low battery alert fuel indicator Lift safety stop if resident or caregiver accidentally gets in the way of the boom arm Audible Alert if the lift detects weight overload Maintenance monitoring: by months or lift cycles Audible Alert that maintenance is due
- Visual battery fuel indicator on Hand Pendant and Lift Control Module
- Auditable low battery alert fuel indicator
- Lift safety stop if resident or caregiver accidentally gets in the way of the boom arm
- Audible Alert if the lift detects weight overload
- Maintenance monitoring: by months or lift cycles
- Audible Alert that maintenance is due
- NEW Six-point cradle with spring clips offer a more secure and wider foot print reducing sling squeeze, shallow height
- offers greater height range
- Can be charged on, or off the lift
- Includes 1 battery and external charger stand

PRODUCTS

Description

FLP500	Levantar Floor lift, power base, 500 lbs
--------	--

SPECIFICATIONS

Description

Battery Pack Weight	7.9 lbs (5.0A), 4.9 lbs (2.9A)
Battery Warranty	1 yr
Front Caster Diameter	3"
Gross Weight	156.6 lbs
Load Capacity Max.	500lbs
Net Weight	132.2 lbs
Overall Height	76.3"
Warranty on Battery	1 yr
Warranty on Electronics	5 yr

ACCESSORIES

Description

FLSCALE	Patient Lift Scale 660lbs capacity, EA
---------	--

ACCESSORIES

Description

13012	Patient Sling,Canvas,1/bx, BX
13019	Chains for Patient Lft,1pr/bx, BX
13023	Patient Lift,New Style,1/cs, CS
13025	Sling,Polystr.for Lift 1/ea, BX
13026	Sling,Polystr.w/Cmde Open1/ea, BX
13057	Replacement Cradle for 13240, 13244 & FL, EA