

drive™



drive™



Bariatric Transfer Board With Hand Holes

Bariatric Transfer Board

Item #	UOM
RTL7047	2/cs

Item #	UOM
RTL7048	2/cs

- For ease of transferring to and from a wheelchair. (Figure A)
- Two cut-out handles provide additional maneuverability.
- Made of sturdy Baltic birch with a clear lacquer finish.
- Ends tapered, corners rounded. Slightly beveled on each end.
- Weight Capacity 600 lbs.
- Size: 35" x 8".

- For ease of transferring to and from a wheelchair. (Figure A)
- Made of sturdy Baltic birch with a clear lacquer finish.
- Ends tapered, corners rounded. Slightly beveled on each end.
- Weight Capacity 600 lbs.
- Size: 29" x 9.5".

specifications

Weight Capacity 600 lbs.

specifications

Weight Capacity 600 lbs.