



Nitro HD Rollator

Item #	Color	UOM
10266HD-BR	Brown	1/ea
<ul style="list-style-type: none">• Attractive, Euro-style design.• Brake cable inside frame for added safety.• Handle height easily adjusts with unique push button.• Back support height easily adjusts with tool-free thumb screw.• Caster fork design enhances turning radius.• Large 10" front casters allow optimal steering and rolling comfort.• Easily folds with one hand to ultra-compact size for storage.• Plastic clip holds the Nitro closed when frame is folded.• Cross-brace design allows for side-to-side folding and added stability.• Seat is durable and comfortable.• Handsome and removable zippered storage bag with unique attachment to keep bag securely in place when Nitro is open or folded.• Lightweight, aluminum frame.		

specifications

Handle Height	36" - 41"
Unit Dimensions	28" (L) x 28.25" (W)
Seat Dimensions	23" (W) x 10" (D)
Seat Height	23.5"
Weight	28.5 lbs.
Weight Limit	450 lbs.
Warranty (Frame)	Limited Lifetime
Warranty (Brake Handle & Housing)	5 Years