



FIGURE A




FIGURE B

- Oxygen tank holder for exclusive use with Drives Nitro rollator
- Accommodates the most commonly used tank sizes (ML6/A tanks, M4/A tanks, M9/C tanks, M6/B tanks, M15/D tanks)
- Conveniently attaches to the Nitros front end without affecting any of the Nitros functionality

**PRODUCTS**

**Description**

10266-40X	Nitro O2 Holder 1/EA,4/CV
10266HD-40X	Nitro HD O2 Holder 1/EA,4/CV

**SPECIFICATIONS**

**Description**

Dimensions 10266-40X	19.5" (W) x 7.5" (L) x 11" (H)
Dimensions 10266HD-40X	24" (W) x 7.5" (L) x 11" (H)
Product Weight	1.7 lbs.