



A



B

Toilet Safety Frame

With Height and Width Adjustable Arms

Item #	Product Description	UOM
12001-4	Assembled	4/cs
12001KD-1	Knocked Down	1/cs

- Anodized aluminum is sturdy and lightweight.
- Plastic armrests provide additional comfort and support.
- Comes with height-adjustable (1" increments) and width-adjustable arms (2 positions). (Figure A)
- Powder-coated aluminum bracket easily attaches frame to toilet.
- Bracket secures safety frame in place. (Figure B)

specifications

Width Between Arms (12001-4)	18.5" (smallest) 22.75" (largest)
Width Between Arms (12001KD-1)	18" (smallest) 25" (largest)
Adjustable Arm Height	25.5" (lowest) 29.5" (highest)
Carton Shipping Weight (12001-4)	18 lbs.
Carton Shipping Weight (12001KD-1)	5 lbs.
Weight Limit	300 lbs.

accessories and parts

12001BR-1	Toilet Safety Frame Bracket	1/ea
-----------	-----------------------------	------