



Toilet Safety Frame

Toilet Safety Rail

Item #	Product Description	UOM
RTL12000	Unassembled	1/cs

- Anodized aluminum is sturdy and lightweight.
- Padded armrests provide additional comfort and support.
- Comes with height-adjustable (1" increments) and width-adjustable arms (2 positions). (Figure A)
- Powder-coated aluminum bracket easily attaches frame to toilet.
- Bracket secures safety frame in place.

specifications

Width Between Arms	19" (smallest) 23.50" (largest)
Adjustable Height	28" (smallest) 32" (highest)
Weight	4 lbs.
Carton Shipping Weight	20 lbs.
Weight Limit	300 lbs.

Item #	UOM
RTL12087	2 retail boxes/cs

- Allows individuals to safely sit or stand.
- Easy, tool-free assembly and installation.
- Comes with soft, padded handles for comfort.
- Width adjusts from 16.5" - 19.5" between handles.

specifications

Width Between Arms	16.5" - 19.5"
Weight	10.5 lbs.
Carton Shipping Weight	24 lbs.
Weight Limit	300 lbs.