



Premium Plastic Raised, Regular/ Elongated Toilet Seat, with Lock

Item #	Product Description	UOM
12014	Without Arms	1/cs

- Designed for individuals who have difficulty sitting down or standing up from the toilet.
- Wide opening in front and back for personal hygiene.
- Lightweight and portable.
- Standard locking mechanism ensures safety. (Figure A)
- Fits most elongated toilets.
- No tools required for installation.
- Easy to clean.

specifications

Seat Dimensions	17" (W) x 16.5" (D) x 4.5" (H)
Weight	4 lbs.
Carton Shipping Weight	6 lbs.
Weight Limit	300 lbs.