

bath safety



Combination Padded Transfer Bench/Commode

Item #	Product Description	UOM
12005KDC-1	Combination Padded Transfer Bench/Commode	1/ea

- Combines transfer bench and commode into one.
- Tool-free assembly of back, legs and arm. (Figure A - C)
- Comfortable cushioned seat and backrest.
- Aluminum A-frame construction is lightweight and provides additional stability.
- Height adjusts in ½” increments with unique “Dual Column” extension legs.
- No exposed hardware to injure patient.
- Extra-large suction cups provide added safety.
- New pinch-free cover allows for push pins to be depressed without pinching fingers. (Figure D)

specifications	
Outside Legs	17.25" (W) x 28.5" (D)
Seat Dimensions	24" (W) x 16" (D)
Seat Height	18.5" - 22.5"
Weight	10 lbs.
Carton Shipping Weight	15 lbs.
Weight Limit	400 lbs.