



## Combination Transfer Bench/ Commode

Tool-free Removable Back and Arm

Item #	Product Description	UOM
12011KDC-2	Knocked Down, Tool-free Removable Back and Arm	2/cs
12011KDC-1	Knocked Down, Tool-free Removable Back and Arm	1/cs

- Tool-free assembly of back, legs and arm. (Figure A - C)
- Combines transfer bench and commode into one.
- Durable blow-molded plastic bench and backrest.
- Aluminum A-frame construction is lightweight and provides additional stability.
- Height adjusts in ½” increments with unique “Dual Column” extension legs.
- No exposed hardware to injure patient.
- Extra-large suction cups provide added safety.
- New pinch-free cover allows for push pins to be depressed without pinching fingers. (Figure D)

### specifications

Outside Legs	31.5" (W) x 22.25" (D)
Seat Dimensions	26.25" (W) x 16.25" (D)
Seat Height	18" (L) x 22" (H)
Weight	11 lbs.
Carton Shipping Weight (-1   -2)	13 lbs.   24 lbs.
Weight Limit	400 lbs.