



Folding Universal Sliding Transfer Bench

Item #	UOM
RTL12075	2/cs

- Transfer Bench slides left and right along frame for easy entry into and out of tub.
- Seat height adjusts in 1/2" increments.
- Folds flat for transport and storage. (Figure A)
- Comes standard with removable soap dish. (Figure B)
- Comes with safety net to catch dropped personal items and prevent bending.
- Comes with seat belt for added safety.
- Retail packaged. (Figure D)

specifications

Seat Dimensions	22.5" (W) x 18.5" (D)
Height	21" - 25"
Weight	16 lbs.
Carton Shipping Weight	38 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime