



Transfer Tub Bench

Knocked down tool-free back, legs and arms

Item #	Product Description	UOM
12011KD-2	Knocked Down, Tool-free Back, Legs, and Arms	2/cs
12011KD-1	Knocked Down, Tool-free Back, Legs, and Arms	1/cs
RTL12011KDR	Knocked Down, Tool-free Back, Legs, and Arms	1/cs

- Tool-free assembly of back, legs and arm. (Figure A - C)
- A-frame construction provides stability.
- Durable blow-molded plastic bench and backrest.
- Height adjusts in ½" increments with unique "Dual Column" extension legs.
- Reversible to accommodate any bathroom.
- Pinch-free lever allows for push pins to be depressed without pinching fingers.

specifications

Outside Legs	32" (W) x 23" (D)
Seat Dimensions	26" (W) x 18.5" - 19.5" (D)
Seat Height	17.5" - 22.5"
Weight	10 lbs.
Carton Shipping Weight (KDR -1 -2)	13 lbs. 13 lbs. 24 lbs.
Weight Limit	400 lbs.