



## Deluxe Folding Walker, Two Button with 5" Wheels

- Easy push-button mechanisms may be operated by fingers, palms or side of hand.
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing.
- Sturdy 1" diameter aluminum construction ensures maximum strength while remaining lightweight.
- Newly designed rear glide cap allows for walker to slide easily and smoothly over most surfaces. (Figure A)
- Vinyl-contoured hand grip.
- Comes with wheels. (Replacement wheels item # 10109).

### specifications

Width	17" inside hand grip   20" inside back legs
Overall Width	24"
Depth (Opened at Base)	17.5"   16.25"
Weight	7.5   7 lbs.
Carton Shipping Weight	30   29 lbs.
Weight Limit	350 lbs.
Warranty	Limited Lifetime