



Flame Red



A



B



C



D

## Deluxe Folding Travel Walker with 5" Wheels and Fold-Up Legs

Item #	Product Description	Color	UOM
RTL10263KDR	Adult	Flame Red Crossbar with Red Frame	1/cs

- Easily folds toll-free to 1/2 the size of a traditional walker (Figure A - B) and can be stored in a "Ready Set Go" Pouch (Figure C)
- Legs are attached with strong elastic cords so legs can be folded. (Figure D - E)
- Comes with 5" front wheels and rear glide skis that are ideal for most surfaces.
- Two side-storage pouches with large pockets for personal items are held in place by hook and loop fastener tabs to both side frames.

specifications	
Width	17" inside hand grip   22" inside back legs
Overall Width	26.5"
Depth (Opened at Base)	19.75"
Height	29.5 - 37"
Weight	8 lbs.
Carton Shipping Weight	10 lbs.
Weight Limit	300 lbs.