



HCPCS
E0141



A

Platform Attachment **Side Walker**

Item #	Product Description	UOM
10219PA	Standard	1/pr
10222PA	Junior	1/pr
10221PA	Pediatric	1/pr

- Height adjustable.
- Can be mounted anywhere along the side bars.
- The handgrips extend, pronate and supinate.
- Can be mounted with an anterior or posterior tilt.
- Brackets can be mounted facing the inside or the outside of the side bars allowing for lateral adjustability.

Item #	Product Description	UOM
10240-1	Adult (Hemi-Walker)	1/ea

- Designed for individuals with the use of only one hand or arm.
- Lighter than a walker and more stable than a cane.
- Folds easily with one hand. (Figure A)
- Ideal for users who do not require a walker but need a wider base for support.

specifications	
Overall Dimensions	18.5" (W) x 16" (D)
Height	28.5" - 35.5"
Weight	4 lbs.
Carton Shipping Weight	5 lbs.
Weight Limit	300 lbs.