



Flame Red



A



B



C



D



Deluxe Folding Travel Walker with 5" Wheels and Fold-Up Legs

Item #	Product Description	Color	UOM
RTL10263KDR	Adult	Flame Red Crossbar with Red Frame	1/cs

- Easily folds tool-free to 1/2 the size of a traditional walker and can be stored in a "Ready Set Go" pouch. (Figure C)
- Legs are attached with strong elastic cords so legs can be folded. (Figure D)
- Comes with 5" front wheels and rear glide skis that are ideal for most surfaces.
- Two side-storage pouches with large pockets for personal items are held in place by Velcro® tabs to both side frames.

specifications

Width	17" inside hand grip 22" inside back legs
Overall Width	26.5"
Depth (Opened at Base)	19.75"
Height	29.5 - 37"
Weight	8 lbs.
Carton Shipping Weight	10 lbs.
Weight Limit	300 lbs.

Steerable Knee Walker

Item #	Product Description	UOM
796	Economy Knee Walker	1/cs

- Knee walker can be steered for increased maneuverability.
- Ideal for individuals recovering from foot surgery, breaks, sprains, amputation and ulcers of the foot.
- New, two-piece padded leg pad provides comfort.
- Deluxe braking system.
- Durable, steel frame with silver vein finish.
- Tool free height adjustment.
- 8" casters are ideal for indoor/outdoor use.
- Limited Lifetime Warranty. Lifetime Warranty on brake cable

specifications

Handle Height	31" - 40"
Unit Dimensions	16" - 30"
Knee Pad Height	16.5" - 20.5"
Pad Dimensions	7" (W) x 15" (D)
Weight	22 lbs.
Carton Shipping Weight	27 lbs.
Weight Limit	300 lbs.