

EasyStand Bantam Sit To Stand Stander, Extra Small	
Base Unit	
Approximate Height Range	28"-40"
Weight Limit	50lbs
Seat Depth Range (from seat pivot to back surface)	7"-12"
Seat Height From Floor (transfer height)	19"
Seat to Footplate Range (from front of seat)	10"-16.5"
Seat to Footplate Range – Inverted Brackets	4.5"-11"
Kneepad Depth Range (from front of seat)	2.5"-8.5"
Kneepad Width	5.5"
Back Height From Seat, Contour & Planar	11"-13"
Back Height From Seat, Contour	12.5"-14.5"
Weight of Basic Unit	54lbs
Frame Footprint – Standard	24.5"x36.5"
Upholstery	Black Nylon/Neoprene
Frame Color	Green, Orange, Purple or White Powder Coat
Lifting Mechanism	Gas Spring Lift
Kneepad (inside access of pad)	2.5"Wx5"H
Footplate Size	7.75"Lx3.25"W
Planar Seat Size	9.5"Wx9.5"D
Planar Back Pad Size	9.5"Wx6.5H
Shadow Tray	
Shadow Tray Size	21"x24"
Shadow Tray Height (from seat)	7"-10"
Weight of Black Molded Shadow Tray Assembly (add 3lbs or 1.4kg for Acrylic)	11.7lbs
Chest Pad Depth Range	0"-11"
Mobile	
Mobile Footprint	24"x38"
Mobile Handrim Height (from seat pivot)	12.75"
Mobile Handrim Height (from floor)	31"
Mobile Handrim Height-Tall (from seat pivot)	16.75"
Mobile Handrim Height-Tall (from floor)	35"
Mobile Wheelbase	31"
Mobile Fore/Aft Handrim Positional Range (from seat pivot)	4.75"-19"
Standard Bantam w/ Mobile Weight	70lbs
Weight of Mobile Option	25lbs
Mobile Drive Chain Tensioning Method	Tension Jack-Screw
Mobile Wheel Locks	1/4 turn rotation to lock/unlock
Mobile Threshold Height	.75" Maximum
Options	
Hip Supports Width Range	6"-9"
Frame Footprint – Swivel Casters	30"x38"
Lifting Mechanism	Manual Hydraulic Lift

Lifting Mechanism	Pow'r Up Lift
Supine Positioning	0-90°
T-Style Handle-Extends (each direction)	7"
Foot Straps – Range Over Arch	2"-9"
Multi-Adjustable Footplates Plantar/Dorsi	+/- 20°
Kneepad Height Range	+/- 1/2"
Multi-Adjustable Footplates Fore/Aft	+/- 3/4"
Multi-Adjustable Footplates Rotation	30° left/right
Multi-Adjustable Kneepads Width Range	5.5"-10"
Hip Support Pad Size	3.5"Wx3.5"H
Hip Supports Height Range	2" (3 positions)
Lateral Support Pad Size	4"Wx6"H
Lateral Supports Width Range	8"-12"
Contoured Seat Size	8.5-9.5"Wx9.5"L
Contoured Back Pad Size	10"Wx6.5"H
Chest Vest	9.5"Lx9W
Head Support Height Range (from seat, small planar only)	17"-24"
Head Support Height Range (all extra small and small contour)	15"-24"
Head Support Depth Range (from backrest)	+1" to -2.75"
Positioning Belt – Hip Circumference	0-23"
Swing-Away Front End Tray Height Range	9.25 to 21.25"
Swing-Away Front End Tray Fore/Aft Range (from seat pivot)	+10 to -11
Locking Caster Diameter	5"