

EasyStand Bantam Sit-To-Stand Stander, Medium	
Base Unit	
Approximate Height Range	48"-66"
Weight Limit	200lbs
Seat Depth Range, Planar (from seat pivot to flat back surface)	15"-20"
Seat Depth Range, Contour (from seat pivot to flat back surface)	15"-20"
Seat (Planar) Height From Floor (transfer height)	22"
Seat (Contour) Height From Floor (transfer height)	25"
Seat to Footplate Range (Planar)	10"-18.5"
Seat to Footplate Range (Planar) -Inverted Brackets	7"-12.5"
Kneepad Depth Range (from seat pivot)	3"-6"
Kneepad Width	8"
Kneepad Dimensions	9"Hx7"W
Back Height From Flat Seat to Top of Back, Short	15"
Back Height From Flat Seat to Top of Back, Tall	19"
Back Height From Contour Seat to Top of Back, Short	16.5"
Back Height From Contour Seat to Top of Back, Tall	20.5"
Weight of Basic Unit	108lbs
Frame Footprint, Standard	24.5"x44.5"
Upholstery	Black Nylon/Neuprene or Molded Urethane
Frame Color	Splatter Midnight
Casters & Wheels (diameter)	5"
Lifting Mechanism	Manual Hydraulic Lift
Seat Pad Size Planar (Width)	15"
Shadow Tray	
Shadow Tray Size	21"x24"
Shadow Tray Height (from top of planar seat)	5.5"-11.25"
Shadow Tray Height (from top of contour seat)	7.25"-12.75"
Shadow Tray Chest Pad Depth Range (from planar back)	2"-10"
Shadow Tray Chest Pad Depth Range (from contour back)	3.5"-11.5"
Weight of Shadow Tray Assembly, Black Molded Tray	10.5lbs
Add for Clear Acrylic Tray	3lbs
Acrylic Tray Size	21"x24"
Oversized Tray Size (Swing-Away Tray Front Only)	22"x29"

Shadow Tray Angle Adjust Range	+40° inclines to -90°
Options	
Frame Footprint (w/ optional Swivel Casters)	28.25"x47"
Lifting Mechanism	Pow'r Up Lift
Supine Positioning	0-90°
T-Style Handle-Extends (each direction)	7"
Adjustable/Removable Actuator Handle Range	+/- 15° from neutral
Foot Straps- fits foot height of	3"-6"
Secure Foot Straps- fits foot height of	3"-6"
Kneepad Height Range	+/- 1.5"
Footplate Size Extra Small	7.75"Lx3.25"W
Footplate Size Small	9.75"Lx4"W
Footplate Size Medium	11.75"Lx5"W
Multi-Adjustable Footplates Plantar/Dorsi	+/- 20°
Multi-Adjustable Footplates Fore/Aft	+/- 0.75" + extra 1.5"
Multi-Adjustable Footplates Rotation	30° left/right
Multi-Adjustable Kneepads Width Range	5"-9"
Elbow Pad, Cutout Width	9.25", 11", or 13"
Elbow Pad, Overall Width	21.5"
Elbow Pad, Extension (back from current table edge)	3.75"
Hip Support Pad Size	4"Hx6.25"W
Hip Supports Width Range	9"-15"
Hip Supports Height Range	+/- 1"
Lateral Support Pad Size	4"Hx6"W
Lateral Supports Width Range	9"-14.5"
Lateral Supports Height Range (tall back)	11"-20"
Lateral Supports Height Range (short back)	10"-17"
Lateral Supports Fore/Aft Adj	3"
Seat Pad Size- Contoured Narrow (Width)	11"-13"
Seat Pad Size- Contoured Wide (Width)	12"-14"
Back Pad Size Planar (Short)	11.5"Hx11.5"W
Back Pad Size Planar (Tall)	15.5"Hx11.5"W
Back Pad Size Contoured (Short)	10.5"Hx11"W
Back Pad Size Contoured (Tall)	14"Hx12"W

Chest Pad Contoured (Large)	9.5"Hx15"W
Chest Strap (Pad Size)	5"Hx8"W
Chest Vest – X-Style (Medium)	9"Lx9.5"W
Chest Vest – X-Style (Large)	11"Lx9.5"W
Chest Vest – Y-Style (Medium)	9"Lx9.5"W
Chest Vest – Y-Style (Large)	11"Lx9.5"W
Head Support Height Range (from seat, short back)	18"-27"
Head Support Height Range (from seat, tall back)	22"-31"
Head Support Depth Range (from backrest)	-2.5" to +1.5"
Head Support-Extra Small	5"Hx8"W
Head Support-Small	6"Hx10"W
Head Support-Large	7"Hx14"W
Positioning Belt- Hip Circumference	to 66"
Positioning Belt with Airline Style Buckle- Hip Circumference	to 66"
Chest Pad Height Range (from seat pivot)	18.5"-33"
Tray Height Range (from seat pivot)	13"-29"
Tray Fore/Aft Range (from seat pivot to standard chest pad)	1"-11"
Frame Footprint (Swing-Away Standard)	27.5"-44"
Frame Footprint (Swing-Away Swivel Front)	29.5"x47.5"