

<b>EasyStand Bantam Sit To Stand Stander, Small</b>	
<b>Base Unit</b>	
Approximate Height Range	36"-54"
Weight Limit	100 lbs
Seat Depth Range (from seat pivot to back surface)	11"-16"
Seat Height From Floor (transfer height)	19"
Seat to Footplate Range (from front of seat)	10"-16.5"
Seat to Footplate Range – Inverted Brackets	4.5"-11"
Kneepad Depth Range (from front of seat)	2.5"-8.5"
Kneepad Width	5.5"
Back Height From Seat, Planar	14"-16"
Weight of Basic Unit	54lbs
Frame Footprint – Standard	24.5"x36.5"
Upholstery	Black Nylon/Neoprene
Frame Color	Green, Orange, Purple or White Powder Coat
Lifting Mechanism	Gas Spring Lift
Kneepad (inside access of pad)	3.5"Wx5"H
Footplate Size	9 3/4"Lx4"W
Planar Seat Size	8"-13"Wx13.5"D
Planar Back Pad Size	8"-12"Wx10"H
<b>Shadow Tray</b>	
Shadow Tray Size	21"x24"
Shadow Tray Height (from seat)	7"-10"
Weight of Black Molded Shadow Tray Assembly (add 3lbs or 1.4kg for Acrylic)	11.7lbs
Chest Pad Depth Range	0"-11"
<b>Mobile</b>	
Mobile Footprint	24"x38"
Mobile Handrim Height (from seat pivot)	12.75"
Mobile Handrim Height (from floor)	31"
Mobile Handrim Height-Tall (from seat pivot)	16.75"
Mobile Handrim Height-Tall (from floor)	35"
Mobile Wheelbase	31"
Mobile Fore/Aft Handrim Positional Range (from seat pivot)	4.75"-19"
Standard Bantam w/ Mobile Weight	70lbs
Weight of Mobile Option	25lbs
Mobile Drive Chain Tensioning Method	Tension Jack-Screw
Mobile Wheel Locks	1/4 turn rotation to lock/unlock
Mobile Threshold Height	.75" Maximum

Options	
Frame Footprint – Swivel Casters	30"x38"
Lifting Mechanism	Manual Hydraulic Lift
Lifting Mechanism	Pow'r Up Lift
Supine Positioning	0-90°
T-Style Handle-Extends (each direction)	7"
Foot Straps – Range Over Arch	2"-9"
Multi-Adjustable Footplates Plantar/Dorsi	+/- 20°
Kneepad Height Range	+/- 1"
Multi-Adjustable Footplates Fore/Aft	+/- 3/4"
Multi-Adjustable Footplates Rotation	30° left/right
Multi-Adjustable Kneepads Width Range	5.5"-10"
Hip Support Pad Size	3.5"Wx3.5"H
Hip Supports Width Range	8"-12"
Hip Supports Height Range	2" (3 positions)
Lateral Support Pad Size	4"Wx6"H
Lateral Supports Width Range	8"-12"
Contoured Seat Size	11-13"Wx14"L
Contoured Back Pad Size	15"Wx8"H
Chest Vest	11"Lx9.5"W
Head Support Height Range (from seat, small planar only)	17"-24"
Head Support Height Range (all extra small and small contour)	15"-24"
Head Support Depth Range (from backrest)	+1" to -2.75"
Positioning Belt – Hip Circumference	0-29"
Swing-Away Front End Tray Height Range	9.25 to 21.25"
Swing-Away Front End Tray Fore/Aft Range (from seat pivot)	+10 to -11
Locking Caster Diameter	5"