

**Bill TO** (Ultimate Medical Supplier Only)

Purchase Order #:	_____
Quote #:	_____
Account #:	_____
Ordered By:	_____
Supplier Company:	_____
Address:	_____
City, State, Zip:	_____
Phone:	_____ Email: _____
Payer Source:	_____
<b>Name of ATP (Required):</b>	_____

**Ship To**

Name:	_____
Address:	_____
City, State, Zip:	_____
Phone:	_____
Tag For:	_____

**Bantam Extra Small User Specifications**

Approximate Height Range	28"-40"
Weight Limit	50 lbs.
Seat Depth Range (seat to back surface)	7"-12"
Seat to Foot Plate Range (from seat pivot)	4.5"-16.5"

**Base** *Note: The DEFAULT items will be added if no selection is made.*

**PT50001-1 EasyStand® Bantam Extra Small Base** \$2,065

**Lift Mechanisms-must select one bolded option below**

- |   |         |
|---|---------|
| <input type="checkbox"/> <b>PT50126 Gas Spring Lift</b> <b>DEFAULT</b>                      | N/C     |
| <input type="checkbox"/> PT50234 Supine (head support required)*                            | \$466   |
| <b>Supine CANNOT be ordered off-unit</b>  |         |
| <input type="checkbox"/> PA5614 Gas Spring Lift Lockout                                     | \$116   |
| <input type="checkbox"/> PT50236 Dual Control (foot and hand control. Push handle required) | \$116   |
| <br>  |         |
| <input type="checkbox"/> <b>PT50020 Manual Hydraulic Actuator with Removable Handle</b>     | \$756   |
| <input type="checkbox"/> PNG50044 Quad Grip Handle Extension                                | \$159   |
| <input type="checkbox"/> PNG50070 T-Style Handle Extension                                  | \$159   |
| <br>  |         |
| <input type="checkbox"/> <b>PT50022 Pow'r Up Lift</b> (battery powered)                     | \$1,265 |

**Frame Color-must select one**

- |   |     |
|---|-----|
| <input type="checkbox"/> Green <b>DEFAULT</b> | N/C |
| <input type="checkbox"/> Purple               | N/C |
| <input type="checkbox"/> White                | N/C |
| <input type="checkbox"/> Orange               | N/C |

**Front Frame-must select one**

- |   |         |
|---|---------|
| <input type="checkbox"/> PT50016-1 Front Frame (5" Front Wheels) <b>DEFAULT</b> | \$163   |
| <input type="checkbox"/> PT50116 Front Frame (3" Front Swivel Casters)          | \$204   |
| <input type="checkbox"/> PT50177 Mobile (not available with supine)             | \$1,588 |

**Trays/Anterior Support-must select one**

- |  |       |
|--|-------|
| <input type="checkbox"/> PA5602 Clear Swing-Away Tray-11"Lx19"W (depth, height & angle adjustable)*    | \$901 |
| <input type="checkbox"/> PA5600 Black Molded Swing-Away Tray-11"Lx19"W <b>DEFAULT</b> *                | \$757 |
| (depth, height & angle adjustable)   |       |
| <input type="checkbox"/> PT50204 Black Molded Tray-11"Lx19"W (depth & height adjustable)*              | \$490 |
| <input type="checkbox"/> PA3038 Hygienic Padded Tray Cover (11"Lx19"W tray only, must select pattern)  | \$213 |
| <input type="checkbox"/> PT50076 Clear Tray-21"Lx24"W (depth, height & angle adjustable)*              | \$892 |
| <input type="checkbox"/> PT50074 Black Molded Tray-21"Lx24"W (depth, height & angle adjustable)*       | \$613 |
| <input type="checkbox"/> PT30093 Hygienic Padded Tray Cover (21"Lx24"W tray only, must select pattern) | \$249 |
| <input type="checkbox"/> PT50246 Mobile Tray-4"Lx12"W  | \$496 |
| (depth and height adjustable, not available w/form to fit backs, lateral supports or elbow stop)       |       |

**Tray Accessories\***

- |   |       |
|---|-------|
| <input type="checkbox"/> PNG50042 Hand Grips (for clear tray only, pair)                    | \$164 |
| <input type="checkbox"/> PT50156 Elbow Pad-9.25" Cutout (21"Lx24"W tray only)               | \$155 |
| <input type="checkbox"/> PT50160 Elbow Pad-11" Cutout (21"Lx24"W tray only)                 | \$155 |
| <input type="checkbox"/> PA5604 Forearm Wings-Black Molded Tray (11"Lx19"W tray only, pair) | \$155 |
| <input type="checkbox"/> PA5606 Forearm Wings-Clear Tray (11"Lx19"W tray only, pair)        | \$155 |

**Knee Supports-must select one**

- |   |       |
|---|-------|
| <input type="checkbox"/> PT50262-1 Slide-In Knee Pads (pair) <b>DEFAULT</b> | \$340 |
| <input type="checkbox"/> PT50264-1 Swing-Away Knee Pads (pair)              | \$563 |

**Knee Pad Size-must select one, inside width**

- |   |     |
|---|-----|
| <input type="checkbox"/> PT50266 2.5"                 | N/C |
| <input type="checkbox"/> PT50268 3.25" <b>DEFAULT</b> | N/C |
| <input type="checkbox"/> PT50270 4.25"                | N/C |

**Knee Pad Accessories**

- |   |       |
|---|-------|
| <input type="checkbox"/> PT50278 Calf Pad | \$147 |
|---|-------|

**Foot Plates-must select one**

- |   |       |
|---|-------|
| <input type="checkbox"/> PT50250-1 Height Adjustable Foot Plates (pair)             | \$275 |
| <input type="checkbox"/> PT50252 Multi-Adjustable Foot Plates (pair) <b>DEFAULT</b> | \$342 |

**Foot Holder Size-must select one**

- |   |     |
|---|-----|
| <input type="checkbox"/> PT50254 6"Lx3"W                      | N/C |
| <input type="checkbox"/> PT50256 7.75"Lx3.25"W <b>DEFAULT</b> | N/C |
| <input type="checkbox"/> PT50258 9.75"Lx4"W                   | N/C |

**Foot Plate Accessories**

- |  |       |
|--|-------|
| <input type="checkbox"/> PT50080 Foot Straps-10"L (length over top of foot, pair)            | \$50  |
| <input type="checkbox"/> PY5636 Foot Straps-15"L (length over top of foot, pair)             | \$50  |
| <input type="checkbox"/> PT50082 Secure Foot Straps-10"L (length over top of foot, two pair) | \$100 |
| <input type="checkbox"/> PY5638 Secure Foot Straps-15"L (length over top of foot, two pair)  | \$100 |

\* not available with mobile

### Seats-must select one

<input type="checkbox"/> PT50042-1 Planar Seat <b>DEFAULT</b>	\$295
<input type="checkbox"/> PT30082 Hygienic Cover for Planar Seat (must select pattern)	\$83
<input type="checkbox"/> PT50044 Form to Fit Seat	\$367
<input type="checkbox"/> PT30083 Hygienic Cover for Form to Fit Seat (must select pattern)	\$83

### Seat Accessories

<input type="checkbox"/> PNG50499 Easy-Adjust Seat Depth- <b>CANNOT be ordered off-unit</b> (roller bearing allows for smoother seat depth adjustment)	\$112
---	-------

### Pelvic Support

<input type="checkbox"/> PT50098 Hip Supports (6"-9" range, pair, not available with form to fit seat)	\$376
<input type="checkbox"/> PA3039 Hygienic Cover for Hip Supports (must select pattern, pair)	\$108
<input type="checkbox"/> PT30066 Positioning Belt (fits hip circumference up to 23", velcro with D-ring)	\$65
<input type="checkbox"/> PT30094 Positioning Belt (fits hip circumference up to 33", velcro with D-ring)	\$65

### Backs-height from seat, must select one

<input type="checkbox"/> PT50050-1 Planar Back-11"-13"H Range <b>DEFAULT</b>	\$295
<input type="checkbox"/> PT30078 Hygienic Cover for PT50050-1 (must select pattern)	\$83
<input type="checkbox"/> PT50238-1 Planar Back-15"-17"H Range	\$295
<input type="checkbox"/> PT30089 Hygienic Cover for PT50238-1 (must select pattern)	\$83
<input type="checkbox"/> PT50052 Form to Fit Back-11"-13"H Range	\$367
<input type="checkbox"/> PT30079 Hygienic Cover for PT50052 (must select pattern)	\$83
<input type="checkbox"/> PT50240 Form to Fit Back-15"-17"H Range	\$367
<input type="checkbox"/> PT30090 Hygienic Cover for PT50240 (must select pattern)	\$83

### Back Accessories

<input type="checkbox"/> PT50212 Push Handle for Supine (required for dual control)	\$159
<input type="checkbox"/> PT50230 Push Handle for Non-Supine (required for dual control)	\$159
<input type="checkbox"/> P80246 Angle Locator	\$40

### Lateral Supports\*

<input type="checkbox"/> PT50008 Lateral Supports-8"-12" Range Curved Pad (for planar back only, pair)	\$376
<input type="checkbox"/> PT50274 Lateral Supports-8"-12" Range Flat Pad (for planar back only, pair)	\$376
<input type="checkbox"/> PT50248 Lateral Supports-5"-9" Range Curved Pad (for planar back only, pair)	\$376
<input type="checkbox"/> PT50276 Lateral Supports-5"-9" Range Flat Pad (for planar back only, pair)	\$376
<input type="checkbox"/> PB3049 Hygienic Cover for Lateral Supports (must select pattern, pair)	\$110
<input type="checkbox"/> PT50216 Lateral Supports with Elbow Stop-8"-12" Range (for planar back only, pair)	\$384
<input type="checkbox"/> PT50218 Lateral Supports with Elbow Stop and Arm Rest-8"-12" Range (for planar back only, pair)	\$425
<input type="checkbox"/> PT50220 Elbow Stop with Arm Rest-8"-12" Range (for planar and form to fit back, pair)	\$384

### Notes:

---



---



---



---

### Upper Body Support

<input type="checkbox"/> PA5672 Padded Positioning Strap	\$159
<input type="checkbox"/> PT50086 X-Style Chest Vest-9.5"Lx9"W	\$297
<input type="checkbox"/> PT50088 X-Style Chest Vest-11"Lx9.5"W	\$297
<input type="checkbox"/> PT50302 Non-Stretch X-Style Chest Vest-9.5"Lx9"W	\$345
<input type="checkbox"/> PNG30318 Non-Stretch X-Style Chest Vest-11"Lx9.5"W	\$345
<input type="checkbox"/> PT50128 High Mount Chest Vest Bracket	\$81

### Head Support-supine requires a head support

<input type="checkbox"/> PT50062 Head Support-5"Hx8"W	\$388
<input type="checkbox"/> PA3035 Hygienic Cover for Head Support-5"Hx8"W (must select pattern)	\$108
<input type="checkbox"/> PT50064 Head Support-6"Hx10"W	\$388
<input type="checkbox"/> PA3036 Hygienic Cover for Head Support-6"Hx10"W (must select pattern)	\$108
<input type="checkbox"/> PT50284 Head Support-Form to Fit-5"Hx10"W	\$388
<input type="checkbox"/> PA3037 Hygienic Cover for Head Support-Form to Fit (must select pattern)	\$108

### Other Accessories

<input type="checkbox"/> P80809 Tool Pouch	\$27
--	------

### Hygienic Cover Patterns (must select one if hygienic covers were ordered)

- Rainbow Unicorn
- Sea Creatures
- Fancy Pants Elephants
- Space Race
- All Smiles Animals

\* not available with mobile

#### (A) Seat Depth Range: 7"-12" (seat to back surface)

Measure the distance between the back of the buttocks and the back of the knee.

#### (B) Seat to Foot Plate Range: 4.5"-16.5" (from seat pivot)

Measure the distance between the foot/shoe and the knee bend.



# AltimateMedical